

# 2000

Multipurpose Patient Chair



800.735.0357  
haag-streit-usa.com



Ultra-durable Upholstery

Swivels 360°

Handrails and Footrest for Safety

Locks In Place with Lever

Corded Footswitch

© 2012 Reliance Medical Products. All Rights Reserved.

# 2000

## Exam Chair for Specialty Procedures

800.735.0357  
haag-streit-usa.com



- **Ergonomic Features**
  - **Handrails and footrest** put patients at ease, and increase safety.
  - **Square base** provides crucial stability during procedures.
  - **Upholstered back** adjusts vertically for patient comfort.
- **Power Base with Corded Footswitch**
  - Corded footswitch raises and lowers the 2000 chair, and can be placed anywhere around the base.
- **Rotation Lock**
  - The versatile 2000 chair swivels 360° and locks in place with a rotation lever.
- **Finishes and Upholstery**
  - Hand-sewn upholstery covers high-density foam that lasts for years.
  - Removable cushions for easy cleaning or replacement.
  - Black upholstery is standard, or choose from 14 standard vinyl colors and endless custom selections.
  - Call Reliance at 1-800-735-0357 for a color card.
- **UL Mark of Quality**
  - The 2000 has passed all rigorous Laboratories tests, achieving CSA, C/US certificate of compliance to assure safety and durability.
- **Tech Specs**
  - Base: 18" square. Seat height: 19¾ lowered, 28¾" raised. Shipping weight: 98 lbs.
  - Electrical: 115 volts | 60 Hz, 2.4 amperes.