

520

Stylish Assisted Recline Tilt Chair



800.735.0357
haag-streit-usa.com



Luxurious Oversized Headrest

Contemporary Design

Upholstered Armrests Rotate Upward

Dual Rotation Locks

Hydraulic Base

© 2012 Reliance Medical Products. All Rights Reserved.

520

Assisted Recline Tilt Exam Chair

800.735.0357
haag-streit-usa.com



● Ergonomic Features

Detached footswitch for hands-free raising/lowering.

Large padded headrest keeps patients comfortable while reclined.

Upholstered armrests rotate up to allow patients to move in and out from either chair side.

Dual-rotation locks allow you to rotate chair 360° from a single position.

Large hinged footrest safely supports patient's full weight, and folds away to save space.

● Effortless Recline to 40°

A manual recline bar, easily accessible from both sides of the chair, controls the assisted tilt of the chair up to 40°.

● Power Hydraulic Base Boosts Efficiency

The 520 chair's corded footswitch allows you to raise and lower the patient while maintaining contact.

● Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning, or replacement.
- Base available in #20 Pearl.
- All materials are fire retardant.
- Choose from 14 standard vinyl colors and endless custom selections – call Reliance at 1-800-735-0357 for a color card.

● CSA Mark of Quality

The 520 has passed all rigorous CSA, C/US tests, affirming its safety and durability and achieving certificate of compliance.

● Tech Specs

Low base height: 21½" lowered, 32½" raised. Base: 22⅞" by 24¼". Depth upright: 43½". Reclined: 59½". Seat: 21½" wide by 15" deep. Armrest distance: 21⅝". Shipping weight: 310 lbs. Standard electrical: 120 volts, 60 Hz, 3.0 amperes. Optional: 230 volts, 50 Hz, 3.0 amperes.