


# The Omega

Omega with sink 103"L x 18-23 1/2"D x 42.5"H



## The Omega

Beautiful and versatile the Omega can be a customized system for the exact needs of your practice. Omega features automatic synchronized room light illumination that responds to your pre-selected levels when you remove the Ophthalmoscope or Retinoscope from wells, lift Binocular Indirect from its cradle, or turn on the Slit Lamp. Reverse any procedure and the room returns to full illumination.

### Standard Features

- ETL approved electronics
- Drug & Rx pad holder
- Four storage drawers
- Built-in tissue dispenser
- Three rechargeable wells
- Unique lazy susan trial lens delivery
- Rheostat-controlled Binocular Indirect hanger/transformer  
(Both are available left or right handed)

# The Performance Leader

Computer Module Version 90" Long x 42.5" High x 18-23.5" Wide



## Omega Optional Computer Module



## Optional Light Box

## Options

- Sink unit
  - Light box drawer
  - Expeditor Call System
  - Computer/scribe area
  - Synchronized fixation, Worth 4 dot and muscle light functions
- Available in panel or mirror set)



## Basic Omega Unit

Omega without sink 74"L x 18-23 1/2"D x 42.5"H